

Easy to see, easy to use FULL COLOR LCD REMOTE CONTROL

High visibility, full color LCD remote control.
Large font sizes with fewer buttons for ease of use.
The vast functions can be easily operated with the up/down and left/right buttons.



Menu button

Various functions can be summoned at the press of a button, including strength adjustment and course changes.

Selection/OK button

Preferred menus and settings can be selected with ease.

Footrest/Reclining/Reset buttons

Adjustments can be made to the reclining and footrest angles.

From the main menu, different courses and massage techniques such as 'Auto Course', 'Basic Knead Massage' and 'Expert Knead Massages' may be selected.

- Select from 10 auto courses.



Recommended Course

·Recommended Course

Whole Body

- Whole Body Rejuvenation Course
- Whole Body Relax Course
- Whole Body Soft Course
- 7 min. (short) course
- Whole Body Stretch Course

Part


- Neck/Shoulder Rejuvenation Course
- Neck/Shoulder Relax Course
- Waist Rejuvenation Course
- Waist Relax Course

- Manual 'Basic Knead Massages' 4 types
- Manual 'Expert Knead Massages' 21 types

**CYBER-RELAX
MASSAGE CHAIR EC-2800**

Specifications

- Measurements (Approx.):
[When the footrest is stored] W710 x D1,340 x H1,150(mm)
[In the reclined position] W710 x D1,890 x H830(mm)
- Minimum width before assembly (Approx.): W710(mm)
- Reclining angle (Approx.): 120~155 degrees
- Weight (Approx.): 72kg, 76kg
- Voltage: AC110~, 220 - 240V~(varies according to the country of use)
- Power consumption: 110W, 130W
- Rated time: 20min.
- Massage chair AS-845 (base domestic Japan model)
Medical Device Licence Number: 224AKBZX00091A01
(Applicable only in Japan)
- *Remote Control Stand
The remote control stand is able to be installed on the left hand side (when seated) of the chair.



Safety
Precautions

- Please read the [User Manual] before use and use this product properly.
- The following persons are strongly recommended not to use this product. Certain physical abnormalities may occur.
 - Persons prohibited from receiving a massage by a medical physician, for example those suffering from thrombosis, severe aneurysm, acute varicose veins, various kinds of dermatitis and skin infections (including inflammations of subcutaneous tissue), etc.
- The following persons must consult a medical physician before using this product as certain physical abnormalities may occur.
 - Persons with malignant tumors. ●Persons who are pregnant or menstruant, persons thought to be pregnant, or those immediately after childbirth.
 - Persons who suffer from paresthesia related to peripheral circulatory disturbance due to diabetes etc. ●Persons with body temperature exceeding 38°C (when strong acute inflammation symptoms, malaise, chills, blood pressure fluctuations are observed, when debilitated) ●Persons with wounds in treatment areas ●Persons who suffer from acute illness (dolorific) ●Persons who suffer from osteoporosis ●Persons who suffer from spinal column fractures ●Persons suffering from inflammation such as sprains and muscle strains ●Persons who suffer from paresthesia ●Persons who suffer from heart trouble ●Persons using a medical electronic device embedded in the body such as a pacemaker susceptible to electromagnetic disturbances
 - Persons with abnormalities in the spine ●Persons who suffer from anemia ●Persons who need to massage a part of the body treated medically in the past
 - Persons who need bed rest ●Persons receiving medical treatment ●Persons who feel an abnormality on the body other than the above
- Please do not allow unattended children or persons unable to declare their own intentions to use this chair. Please do not allow infants near the chair.

●Concerning the use of this catalogue

- The material contained in this catalogue are based on domestic Japanese Pharmaceutical Affairs Laws. The claims made in this catalogue may not be valid or relevant in your country.
- Fujiiryoki will not be held responsible for the use of this catalogue outside of Japan, and all responsibility will lie with the user in such circumstances.
- Product images shown may differ from actual product appearance.
- Designs and specifications are subject to change without notice.

FUJIIRYŌKI

Fuji Medical Instruments Mfg. Co., Ltd.
Overseas Sales Department


14th Floor, Oe Building, 1-22, Noninbashi 1-chome, Chuo-ku, Osaka, 540-0011, Japan
Tel: +81-6-4793-0745 Fax: +81-6-4793-0781

FUJIIRYŌKI

MASSAGE CHAIR EC-2800

HYBRID

**Large
Knead Balls**



**Digital Mecha
Control**



CYBER-RELAX

THE PIONEER IN MASSAGE CHAIRS

Since 1954

HYBRID

A new massage experience!!

Introducing Fujiiryoki's new hybrid mecha knead ball system.

Large 9cm knead balls with state of the art digital mecha control.

By re-enacting a 'palm' massage (where the base of the palm is used to give a massage), we have been able to realize a gentle, gradual massage over a wider, larger area.

Covering the shoulders to the waist, with just the right amount of pressure, the 9cm knead balls knead and loosen stiff muscles with minimal discomfort for a experience only available from Fujiiryoki's new hybrid mecha knead ball system.



Large,
wide
profile

EC-2800

Original size
(EC-2700)

The larger knead balls stimulate a wider area compared with traditional knead balls, re-enacting a palm massage. *Comparison based on Fujiiryoki products.

Large sized knead balls

Diameter
9cm (approx.)
×
Thickness
2.5cm (approx.)



Fujiiryoki's original digital mecha control

Independent
strength
adjustability
for each part
of the body



SHOULDER

Strength
Adjustment
8 Levels



BACK

Strength
Adjustment
2 Levels



WAIST

Strength
Adjustment
8 Levels



Expert
Knead
Massages

25 different massage techniques
(including basic knead massages)
designed for different purposes and
areas of the body.

Discover the distinctive characteristics of the large knead balls around the shoulders and back and especially the wide area massages around the waist.



LEGS Large airbags provide a 'wrapping' massage for the calves, the back of the calves, ankles, foot soles, and the area under the knees.

Airbags Strength Adjustment **2 Levels**

Leg Air Massage



12 airbags wrap the entire legs, from the calves to the feet soles, applying repeated pressure followed by a releasing action to promote blood circulation.

Foot Stretch Function



The airbags firmly hold the legs and pull down to create a stretching massage.

Spring Slide Function



The lower section of the 2-step footrest can be extended to a maximum 16cm to allow adjustment for different sized people.



SHOULDER~ARMS

No matter the body shape or size, a comfortable massage for the shoulders and upper arms, from the arms to the fingertips.

Knead Balls Strength Adjustment **8 Levels**

Airbags Strength Adjustment **2 Levels**

Shoulder Muscle Care



Focus on kneading and loosening up of the muscles, moving up from the shoulders to the base of the neck. Particularly effective in relieving stiffness and pain of the muscles.

Shoulder Care Massage



Airbags rising up from the back gently wrap the shoulders from behind, promoting blood circulation around the shoulder area.

Arm Care Massage



4 airbags, located above and to the side, stimulate points from the arms to the fingertips, relieving stiffness and fatigue of the arms.



Knead Ball Position Adjustment Function

Before commencing the massage, the knead balls will automatically detect the position of the shoulder/neck and make adjustments from 9 different positions to provide the perfect position for different people of different shapes and sizes.




BACK

Effectively kneads and loosens the area around the back and the shoulder blades to relieve stiffness and fatigue.

Knead Balls Strength Adjustment **2 Levels**

NEW


Spine Stretch



A massage which stretches the muscles in the back with side movements moving from the center of the back to the outer areas. Kneads and loosens to relieve pain from stiffness.


NEW

Spine Care



The knead balls travel up and down and stretch the muscles in the entire back. Especially recommended for persons who are troubled by stiffness, pain, and fatigue in their back area caused by work or household labor.

Shoulder Blade Line



The knead balls trace the shoulder blade line and gently massage the stiff muscles behind the shoulder blades caused by desk work and spending long hours in the same posture.



Strength Adjustment 2 Levels

Independent Air Strength Area Adjustability

The strength of the air massage is adjustable to 2 levels on each part of the body.



WAIST

Firmly kneads the waist and hip areas, with the airbags wrapping and gently massaging the pelvic area.


Knead Balls Strength Adjustment **8 Levels**

Airbags Strength Adjustment **2 Levels**

KIWAMI-Knead, KIWAMI-Tapping



Recommended for those who suffer from nerve or muscle pain, or general fatigue in the waist/hip areas. Effective pinpoint kneading and tapping for the waist.



Pelvic Air Massage

The combination of airbags to the side of the thighs and in the seat wrap and gently massage the muscles around the pelvic area from the sides and from below. Staggered movements of the airbags also provide a stretching action for the area around the waist.



* Other kneading and tapping massages for a total 25 different types of basic kneading techniques.

*Images for illustration purposes only.